

# AI/ML-Based Personalized Diet Plan Generator from Medical Reports

## Milestone 4: Diet Plan Generation and UI Integration

### Milestone Objectives

The objective of Milestone 4 was to integrate all previously developed components into a **fully functional end-to-end system** capable of generating personalized diet plans based on both structured and unstructured medical inputs.

As per the project roadmap, this milestone aimed to:

- Combine ML-based health analysis and AI/NLP-based text interpretation
- Generate personalized diet plans using integrated insights
- Implement a user-friendly frontend interface using React
- Enable export of diet plans in PDF and JSON formats
- Deliver a fully functional prototype tested with sample medical reports

### System Integration Overview

By Milestone 4, the system combined three major components:

1. **Structured Data Analysis (Milestone 2)**
  - ML model (XGBoost) predicts diabetes risk
  - Provides probability score and risk level
  - Uses ROC-AUC as reliability metric
2. **Unstructured Text Analysis (Milestone 3)**
  - OCR extracts text from uploaded prescriptions
  - NLP interprets dietary instructions
  - Extracted rules mapped to diet guidelines
3. **Diet Plan Generation Engine (Milestone 4)**
  - Combines numeric risk analysis + extracted dietary rules
  - Generates structured daily meal plans
  - Displays results in a modern frontend UI

This milestone completes the full pipeline from **input** → **analysis** → **interpretation** → **personalized output**.

## **Personalized Diet Plan Generation Logic**

The diet generation module uses:

- ML-based diabetes risk level
- Extracted prescription instructions
- Predefined nutritional rules
- Balanced macronutrient distribution

Based on the combined insights, the system generates a complete daily meal plan structured into:

- **Breakfast**
- **Lunch**
- **Dinner**
- **Recommended Snacks**

Each meal includes:

- Food item name
- Description
- Calorie value
- Macronutrient indicators (carbohydrates, protein, fat)

The plan adapts based on:

- Risk level (Low / Moderate / High)
- Dietary restrictions (e.g., avoid sugar, low carb, high fiber)
- General health recommendations

## **Frontend Implementation (React UI)**

The frontend was implemented using React to provide an interactive and user-friendly interface.

**Core UI Features:**

- **Input form for health parameters**
- **Prescription image upload**
- **Real-time risk prediction display**
- **Dynamic meal plan rendering**
- **Color-coded risk visualization**

- **Export options**

The UI ensures that users can easily:

- **Enter medical data**
- **Upload prescriptions**
- **View personalized recommendations**
- **Download results**

### **Export Functionality**

To improve usability and portability, the system supports:

#### **1. PDF Export**

- Allows users to download their personalized meal plan
- Useful for offline access or medical consultation

#### **2. JSON Export**

- Structured output format
- Enables future integration with:
  - Mobile apps
  - Nutrition tracking tools
  - Healthcare systems

### **End-to-End Prototype Validation**

The system was tested using:

- Sample numeric medical inputs
- Sample prescription images
- Extracted dietary instructions

#### **Validation Results:**

- Risk prediction successfully generated
- OCR extracted readable text
- NLP mapped dietary instructions correctly
- Personalized meal plan generated dynamically
- Export functionality worked successfully

The system achieved a **fully functional end-to-end prototype** meeting the milestone criteria.

**System Architecture Completion**

By the end of Milestone 4, the following modules are fully integrated:

- Medical Data Input Module
- OCR Extraction Module
- NLP Interpretation Engine
- ML-Based Health Risk Analysis
- Personalized Diet Generation Engine
- React-Based Frontend Interface
- Export (PDF & JSON) Module

This completes the core architecture of the AI Diet Planner system.

**UI Demonstration**

**Figure: End-to-End AI Diet Planner Interface with Generated Meal Plan**

The screenshot demonstrates:

- Extracted medical information
- Generated breakfast, lunch, dinner, and snack recommendations
- Calorie and nutrient indicators
- Export options

**OPTION 1: Manual Entry**

AI Diet Planner

AI-Powered Diabetes Risk & Personalized Diet Recommendation System

OPTION 1: MANUAL ENTRY

Patient Health Parameters

Enter health metrics manually if you know the values

Pregnancies

0-17

Glucose Level (mg/dL)

0-1000

Blood Pressure (mmHg)

0-300

Skin Thickness (mm)

0-100

Insulin Level (µIU/mL)

0-900

BMI (kg/m²)

0-100

Diabetes Pedigree Function

0-3

Age (years)

1-150

PREDICT RISK

GENERATE FULL MEAL PLAN

CLEAR FORM

OPTION 1: MANUAL ENTRY

Patient Health Parameters

Enter health metrics manually if you know the values

|             |                            |                       |                     |                        |
|-------------|----------------------------|-----------------------|---------------------|------------------------|
| Pregnancies | Glucose Level (mg/dL)      | Blood Pressure (mmHg) | Skin Thickness (mm) | Insulin Level (μIU/mL) |
| 1           | 400                        | 100                   | 30                  | 400                    |
| BMI (kg/m²) | Diabetes Pedigree Function | Age (years)           |                     |                        |
| 40          | 1                          | 45                    |                     |                        |

PREDICT RISK

GENERATE FULL MEAL PLAN

CLEAR FORM

Prediction Results

RISK LEVEL:

High Risk

DIABETES PROBABILITY:

78.7%

MODEL ROC-AUC SCORE:

0.81

RECOMMENDATION:

Follow a low-GI diet and consult a doctor

Personalized Meal Plan

Based on your health parameters and diabetes risk

Risk Level:

High Risk

Diabetes Probability:

78.7%

MONDAY

BREAKFAST Oatmeal with berries and nuts 350 cal

LUNCH Grilled chicken salad 420 cal

SNACK Whole grain crackers with cheese 140 cal

DINNER Baked salmon with steamed vegetables 450 cal

DAILY TOTAL ~1360 cal

TUESDAY

BREAKFAST Vegetable omelet with whole wheat toast 320 cal

LUNCH Quinoa bowl with roasted vegetables 380 cal

SNACK Greek yogurt 120 cal

DINNER Grilled chicken with brown rice 480 cal

DAILY TOTAL ~1300 cal

WEDNESDAY

BREAKFAST Greek yogurt with chia seeds 280 cal

LUNCH Lentil soup with vegetables 350 cal

SNACK Fresh fruit 80 cal

DINNER Turkey meatballs with zucchini noodles 420 cal

DAILY TOTAL ~1130 cal

THURSDAY

BREAKFAST Oatmeal with berries and nuts 350 cal

LUNCH Grilled chicken salad 420 cal

SNACK Handful of nuts 160 cal

DINNER Baked salmon with steamed vegetables 450 cal

DAILY TOTAL ~1380 cal

FRIDAY

BREAKFAST Vegetable omelet with whole wheat toast 320 cal

LUNCH Quinoa bowl with roasted vegetables 380 cal

SNACK Vegetable sticks with hummus 100 cal

DINNER Grilled chicken with brown rice 480 cal

DAILY TOTAL ~1280 cal

SATURDAY

BREAKFAST Greek yogurt with chia seeds 280 cal

LUNCH Lentil soup with vegetables 350 cal

SNACK Whole grain crackers with cheese 140 cal

DINNER Turkey meatballs with zucchini noodles 420 cal

DAILY TOTAL ~1190 cal

SUNDAY

BREAKFAST Oatmeal with berries and nuts 350 cal

LUNCH Grilled chicken salad 420 cal

SNACK Whole grain crackers with cheese 140 cal

DINNER Baked salmon with steamed vegetables 450 cal

DAILY TOTAL ~1360 cal

EXPORT AS PDF

EXPORT AS JSON

## OPTION 2: Upload Doctor Prescription

### OPTION 2: UPLOAD REPORT

## Upload Doctor Prescription

Upload a medical report to automatically extract health parameters

Selected: MUM-0425-PA-0004291\_LAB REPORT\_26-04-2025\_0523-21\_PM@G.pdf\_page\_1.png

 ANALYZE & EXTRACT

 CLEAR FILE

Selected: MUM-0425-PA-0004291\_LAB REPORT\_26-04-2025\_0523-21\_PM@G.pdf\_page\_1.png (624.09 KB)

## Extraction Results

### Extracted Parameters

Blood Pressure:

111

### Extracted Text

RENAL PANEL - I  
Specimen: Serum/ Plasma  
BUN (Urease/GLDH) 11.60 mg/dl [8.00-23.00]  
SERUM CREATININE (modJaffe) 0.64 mg/dl [0.50-0.90]  
\*eGFR(Calculated) 94.6 ml/min/1.73sq.m [>60.0]  
SERUM SODIUM (Indirect ISE) 139.5 mmol/l [136.0-145.0]  
SERUM POTASSIUM (Indirect ISE) 4.06 mmol/l [3.50-5.10]  
SERUM CHLORIDE (Indirect ISE) 102.5 mmol/L [98.0-107.0]  
SERUM BICARBONATE(Enzymatic PEP-MD) 23.4 mmol/L [21.0-31.0]  
UREA 24.8 mg/dl [17.0-43.0]

### Diet Rules

- ✓ Soft, easy-to-digest foods
- ✓ High protein for recovery and healing
- ✓ Nutrient-dense foods to boost immunity
- ✓ Focus on: Zinc, Vitamins A, C, D, Protein
- ✓ Avoid: Processed junk food



## 7-Day Personalized Meal Plan

Complete weekly meal plan with breakfast, lunch, snack, and dinner

### MONDAY

#### BREAKFAST

Oatmeal with berries & almonds

350 cal 12g protein

#### LUNCH

Grilled chicken breast with brown rice

520 cal 45g protein

#### SNACK

Apple with almonds

180 cal 6g protein

#### DINNER

Baked lean fish with sweet potato

550 cal 48g protein

Daily Total: 1600 cal

### TUESDAY

#### BREAKFAST

Scrambled eggs with whole wheat toast

300 cal 18g protein

#### LUNCH

Salmon fillet with steamed broccoli

480 cal 42g protein

#### SNACK

String cheese & crackers

150 cal 8g protein

#### DINNER

Grilled chicken with quinoa

580 cal 50g protein

Daily Total: 1510 cal

### WEDNESDAY

#### BREAKFAST

Greek yogurt with granola

280 cal 15g protein

#### LUNCH

Turkey sandwich with lettuce

420 cal 35g protein

#### SNACK

Carrot sticks with hummus

120 cal 4g protein

#### DINNER

Beef stir-fry with vegetables

520 cal 46g protein

Daily Total: 1340 cal

### THURSDAY

#### BREAKFAST

Smoothie bowl with banana & chia

320 cal 10g protein

#### LUNCH

Chickpea Buddha bowl

480 cal 18g protein

#### SNACK

Mixed nuts

170 cal 6g protein

#### DINNER

Vegetable curry with lentils

480 cal 20g protein

Daily Total: 1450 cal

### FRIDAY

#### BREAKFAST

Whole wheat pancakes with honey

340 cal 14g protein

#### LUNCH

Tuna salad with olive oil

450 cal 38g protein

#### SNACK

Protein bar

200 cal 10g protein

#### DINNER

Turkey meatballs with pasta

520 cal 42g protein

Daily Total: 1510 cal

### SATURDAY

#### BREAKFAST

Vegetable omelette with mushrooms

280 cal 16g protein

#### LUNCH

Tofu stir-fry with vegetables

380 cal 28g protein

#### SNACK

Berries with yogurt


140 cal 7g protein

#### DINNER

Baked chicken breast with veggies

500 cal 52g protein

Daily Total: 1300 cal



## Key Learnings

Through Milestone 4, the following skills and concepts were strengthened:

- Full-stack AI application development
- Integration of ML + NLP pipelines
- Backend-frontend API communication
- User-centered UI design
- Deployment-ready architecture design
- End-to-end healthcare system prototyping



## **Conclusion**

Milestone 4 successfully integrated machine learning-based health analysis and NLP-based prescription interpretation into a complete, functional system capable of generating personalized diet plans. The frontend implementation provides an intuitive interface, while the backend ensures intelligent processing and reliable results.

The system now operates as a **fully functional end-to-end prototype**, tested with sample medical data and prescription reports. This milestone marks the completion of the core implementation phase of the AI Diet Planner project.

**SUBMITTED BY:**

**Sabiha Anjum**